

Jeio Tech 55 Litre Shaking Water Bath (BS-31)

Product Code: 44312

Brand: Jeio Tech

The Jeio Tech 55 Litre Shaking Water Baths (BS-31) is a reciprocal shaking and heating bath that provides precise & constant temperature control and is ideal for sample incubation conditions for DNA, tissue, or cell culture applications.

This water bath is supplied with a hinged gable cover and spring wire rack.

Features:

Performance

- Temperature range from ambient +5°C to 100°C.
- Reciprocating shaking up to 180rpm while maintaining uniform temperature.
- Microprocessor PID control/temperature calibration/automatic tuning/automatic run.
- Separate control keys for shaker mechanism allows the unit also to be used as a heating bath.

Convenience

- Multiple LED displays (0.1°C and 1rpm resolution) for various values.
- Dual wait on/off timer modes (1 min. to 99 hr. 59 min.).
- Minimal shock to the samples thanks to the shaker's gentle feature.
- Spring wire rack is easy to equip and the height of the rack is easily adjustable.
- Silent and anti-vibration BLDC motor.
- SUS 304 spring wire rack (two tier rack with flexible crossed spring) provides various ways of glassware settings.

Safety

- Isolated design of agitator, heater, and sensor by a baffle plate to minimise the turbulence in the bath fluid and also to protect users and samples against direct contact.
- Complete safety protection system with warning alarms.



Specifications:

Bath Volume (litres)		55
Working Temperature Range (°C)		Amb. +5 to 100
Temperature Stability (°C)		± 0.2
Shaking System	Motion/Speed Range (rpm)	Reciprocating/20 to 180
	Recip. Stroke Length (mm)	30 (25 and 35 available)
Dimensions	Bath Opening/Depth (w x l, d)	629 x 350, 235mm
	External (w x d x h)	860 x 490 x 355mm
	Net Weight (kg)	42
Electrical Requirements (230V, 50Hz)		8.9A

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com