

Deltalab 250ml Wide Neck Containers with Star-Shaped Caps

Product Code: 44440000

Brand: Deltalab

Made of high density polyethylene, Deltalab's 250ml Wide Neck Containers with Star-Shaped Caps have tagging points on the neck and cap for a tamper proof closure or for labelling. These laboratory containers are supplied uncapped and are made with materials suitable for alimentary use.



Specifications

Volume (ml)	250
Internal Mouth (Ø mm)	37
Body (Ø mm)	65
Height Without Lid (mm)	128
Body Weight (g)	24
Quantity	120

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com