

Deltalab 1000ml Wide Neck Containers with Star-Shaped Caps

Product Code: 44460000

Brand: Deltalab

Made of high density polyethylene, Deltalab's 1000ml Wide Neck Containers with Star-Shaped Caps have tagging points on the neck and cap for a tamper proof closure or for labelling. These laboratory containers are supplied uncapped and are made with materials suitable for alimentary use.



Specifications

Volume (ml)	1000
Internal Mouth (Ø mm)	59.4
Body (Ø mm)	102
Height Without Lid (mm)	183.2
Body Weight (g)	90
Quantity	55

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com