

Deltalab 60ml White Security Screw Cap Containers with Blue Caps

Product Code: 444602801

Brand: Deltalab

Deltalab's 60ml White Security Screw Cap Containers with Blue Caps are high density polyethylene containers and are ideal for liquids.

These laboratory containers have a screw cap and provide leakproof closure with their child-proof system and split seal.

The cap is compound by two pieces that turn around themselves avoiding the cap to be opened, working as a security closure avoiding the accidental manipulation. To open the cap, push it down and turn it. The sealed caps need to be unscrewed to break the seal, showing that it has not been manipulated.



Specifications

Capacity (ml)	60
Cap Colour	Blue
Body Colour	White
Height (mm)	87.5
Base (Ø mm)	37.6
Body Weight (g)	9.5
Quantity	600

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com