

Deltalab 60ml Natural Security Screw Cap Containers with White Caps

Product Code: 444603300

Brand: Deltalab

Deltalab's 60ml Natural Security Screw Cap Containers with White Caps are high density polyethylene containers and are ideal for liquids.

These laboratory containers have a screw cap and provide leakproof closure with their child-proof system and split seal.

The cap is compound by two pieces that turn around themselves avoiding the cap to be opened, working as a security closure avoiding the accidental manipulation. To open the cap, push it down and turn it. The sealed caps need to be unscrewed to break the seal, showing that it has not been manipulated.



Specifications

Capacity (ml)	60
Cap Colour	White
Body Colour	Natural
Height (mm)	93.7
Base (Ø mm)	37.6
Body Weight (g)	9.5
Quantity	600

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com