

## Deltalab 25ml Square Containers with Star-Shaped Caps

**Product Code:** 444608

**Brand:** Deltalab

Deltalab's 25ml Square Containers with Star-Shaped Caps have tagging points on the neck and cap for a tamper proof closure or for labelling. The containers are made of high density polyethylene and the star-shaped screw caps are made of polypropylene.

These labware containers have moulded-in graduations and are made with materials suitable for alimentary use.

The containers are supplied plugged and capped. They are not recommended for liquid transport.



### Specifications

Capacity (ml)	25
Graduation Each (ml)	12.5
Dimensions (mm)	32 x 34 x 52
Neck (Ø mm)	18
Body Weight (g)	6
Quantity	500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)