

## Deltalab 50ml Square Containers with Star-Shaped Caps

**Product Code:** 444609

**Brand:** Deltalab

Deltalab's 50ml Square Containers with Star-Shaped Caps have tagging points on the neck and cap for a tamper proof closure or for labelling. The containers are made of high density polyethylene and the star-shaped screw caps are made of polypropylene.

These labware containers have moulded-in graduations and are made with materials suitable for alimentary use.

The containers are supplied plugged and capped. They are not recommended for liquid transport.

### Specifications

<b>Capacity (ml)</b>	50
<b>Graduation Each (ml)</b>	10
<b>Dimensions (mm)</b>	38 x 38 x 70
<b>Neck (Ø mm)</b>	23
<b>Body Weight (g)</b>	15
<b>Quantity</b>	400



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)