

Deltalab 250ml Square Containers with Star-Shaped Caps

Product Code: 444611

Brand: Deltalab

Deltalab's 250ml Square Containers with Star-Shaped Caps have tagging points on the neck and cap for a tamper proof closure or for labelling. The containers are made of high density polyethylene and the star-shaped screw caps are made of polypropylene.

These labware containers have moulded-in graduations and are made with materials suitable for alimentary use.

The containers are supplied plugged and capped. They are not recommended for liquid transport.

Specifications

Capacity (ml)	250
Graduation Each (ml)	50
Dimensions (mm)	57 x 60 x 110
Neck (Ø mm)	33
Body Weight (g)	32
Quantity	300



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com