

Deltalab 1000ml Square Containers with Star-Shaped Caps

Product Code: 444613

Brand: Deltalab

Deltalab's 1000ml Square Containers with Star-Shaped Caps have tagging points on the neck and cap for a tamper proof closure or for labelling. The containers are made of high density polyethylene and the star-shaped screw caps are made of polypropylene.

These labware containers have moulded-in graduations and are made with materials suitable for alimentary use.

The containers are supplied plugged and capped. They are not recommended for liquid transport.

Specifications

| | |
|-----------------------------|----------------|
| Capacity (ml) | 1000 |
| Graduation Each (ml) | 100 |
| Dimensions (mm) | 80 x 103 x 176 |
| Neck (Ø mm) | 58 |
| Body Weight (g) | 97 |
| Quantity | 80 |



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com