

Deltalab 2000ml Square Containers with Star-Shaped Caps

Product Code: 444614

Brand: Deltalab

Deltalab's 2000ml Square Containers with Star-Shaped Caps have tagging points on the neck and cap for a tamper proof closure or for labelling. The containers are made of high density polyethylene and the star-shaped screw caps are made of polypropylene.

These labware containers have moulded-in graduations and are made with materials suitable for alimentary use.

The containers are supplied plugged and capped. They are not recommended for liquid transport.

Specifications

Capacity (ml)	2000
Graduation Each (ml)	100
Dimensions (mm)	100 x 134 x 210
Neck (Ø mm)	58
Body Weight (g)	177
Quantity	40



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com