

## Deltalab 2000ml Wide Neck Containers with Star-Shaped Caps

**Product Code:** 44470000

**Brand:** Deltalab

Made of high density polyethylene, Deltalab's 2000ml Wide Neck Containers with Star-Shaped Caps have tagging points on the neck and cap for a tamper proof closure or for labelling. These laboratory containers are supplied uncapped and are made with materials suitable for alimentary use.



### Specifications

<b>Volume (ml)</b>	2000
<b>Internal Mouth (Ø mm)</b>	59.4
<b>Body (Ø mm)</b>	123.5
<b>Height Without Lid (mm)</b>	244.5
<b>Body Weight (g)</b>	128
<b>Quantity</b>	30

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)