

Deltalab Yellow Tissue Embedding Cassettes for Biopsy (Histosette II) - Lid Pre-Mounted

Product Code: 445143

Brand: Deltalab

Made of acetal polymer, Deltalab's Yellow Tissue Embedding Cassettes for Biopsy (Histosette II) - Lid Pre-Mounted have 224 holes and 1mm square openings to maximize fluid exchange and to ensure proper drainage.


The lids are available pre-mounted.

Features

- Hinged lid can be opened and closed as often as necessary
- Supplied assembled, saving time and making storage easier
- Air vents allowing more efficient filling with paraffin
- Easy handling with one hand
- Safety closure
- Writing area slanted at a 45° angle, suitable to be used with most labelling instruments

Specifications

Colour	Yellow
External Dimensions (mm)	41 x 28 x 6
Quantity	1500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com