

## Deltalab Yellow Tissue Embedding Cassettes (Histosette I)

**Product Code:** 454143

**Brand:** Deltalab


Deltalab's Yellow Tissue Embedding Cassettes (Histosette I) have a writing area at 30° and 64 cuts in the base and lid.

### Features

- Made of polyacetal
- Air vents allowing more efficient filling with paraffin
- Safety closure
- With one-piece integral lid
- To lock the cassette, just snap apart the lid, turn it round and lock it into the base of the cassette
- They may be opened and closed as many times as required

### Specifications

<b>Colour</b>	Yellow
<b>External Dimensions (mm)</b>	41 x 28 x 6
<b>Quantity</b>	1500



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)