

Deltalab Green Tissue Embedding Cassettes (Histosette I)

Product Code: 454147

Brand: Deltalab


Deltalab's Green Tissue Embedding Cassettes (Histosette I) have a writing area at 30° and 64 cuts in the base and lid.

Features

- Made of polyacetal
- Air vents allowing more efficient filling with paraffin
- Safety closure
- With one-piece integral lid
- To lock the cassette, just snap apart the lid, turn it round and lock it into the base of the cassette
- They may be opened and closed as many times as required

Specifications

Colour	Green
External Dimensions (mm)	41 x 28 x 6
Quantity	1500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com