

Deltalab Round Jugs with Handle (0.5 Litres)

Product Code: 470000

Brand: Deltalab

Made of polypropylene, Deltalab's Round Jugs with Handle (0.5 Litres) are autoclavable up to 121°C and are made with materials suitable for alimentary use.

They feature moulded-in graduations.



Specifications

Capacity (litres)	0.5
Graduation (ml)	25
Diameter (mm)	80
Height (mm)	120
Material	Polypropylene
Case Quantity	60

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com