

## Deltalab 5 Litre Industrial Jugs

**Product Code:** 470007

**Brand:** Deltalab

Deltalab's 5 Litre Industrial Jugs are made of thick rigid polyethylene and have a natural colour. They are recommended for industrial use and made with materials suitable for alimentary use.



## Specifications

<b>Capacity (L)</b>	5
<b>Mouth (Ø mm)</b>	190
<b>Height (mm)</b>	240
<b>Weight (ut. g)</b>	400
<b>Quantity</b>	1

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)