


Deltalab Measuring Cylinders (1000ml)

Product Code: 481000

Brand: Deltalab

Deltalab's Measuring Cylinders (1000ml) have stable pentagonal bases, along with printed graduations for clear volume evaluation. Each cylinder is autoclavable up to 121°C and in compliance with ISO 6706 and BS 5404 standards.



Specifications:

Capacity (ml)	1000
Graduation (ml)	200
Subdivision (ml)	20
Overall Diameter (mm)	80
Height (mm)	530
Material	Polypropylene
Sterility	Non-Sterile
Quantity	6

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com