


## Deltalab Measuring Cylinders (2000ml)

**Product Code:** 482000

**Brand:** Deltalab

Deltalab's Measuring Cylinders (2000ml) have stable pentagonal bases, along with printed graduations for clear volume evaluation. Each cylinder is autoclavable up to 121°C and in compliance with ISO 6706 and BS 5404 standards.



### Specifications:

<b>Capacity (ml)</b>	2000
<b>Graduation (ml)</b>	200
<b>Subdivision (ml)</b>	20
<b>Overall Diameter (mm)</b>	80
<b>Height (mm)</b>	530
<b>Material</b>	Polypropylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	6

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)