

Deltalab Graduated Three Spouts Beakers (50ml)

Product Code: 490050

Brand: Deltalab

Made of ultra clear polypropylene, Deltalab's Graduated Three Spouts Beakers (50ml) are autoclavable up to 121°C. These disposable beakers feature moulded-in graduations.

Specifications:

Volume (ml)	50
Minimum Order Quantity	10
Material	Polypropylene
Case Quantity	100



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com