


## Deltalab Graduated Beakers (2000ml)

**Product Code:** 492000

**Brand:** Deltalab

Ideal for general laboratory use, Deltalab's Graduated Beakers (2000ml) are made of durable polypropylene for safe use. These lab beakers are autoclavable up to 121°C, have clearly printed graduations, and are in compliance with ISO 6706 and BS 5404 standards.



### Specifications

<b>Capacity (ml)</b>	2000
<b>Diameter (mm)</b>	135
<b>Height (mm)</b>	184
<b>Material</b>	Polypropylene
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	4

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)