

Deltalab 250ml Screw Cap Containers with Seal and Dropper

Product Code: 53250

Brand: Deltalab

Made of white polypropylene, Deltalab's 250ml Screw Cap Containers with Seal and Dropper have a watertight screw closing system with dropper plug and blue security seal.

They have a hermetical closure and are supplied uncapped.

Specifications

Volume (ml)	250
Base (Ø mm)	50.8
Height (mm)	167.3
Body Weight (g)	21.50
Mouth Diameter (mm)	10.2
Quantity	200



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com