

## Deltalab Pink Tissue Embedding Cassettes

**Product Code:** 554140


**Brand:** Deltalab

Deltalab's Pink Tissue Embedding Cassettes are made of acetal polymer and have integral lids.

Each cassette has a small tab on the lid for easy opening and closing.

### Features

- Air vents allowing more efficient filling with paraffin
- Safety closure
- With one-piece integral lid
- To lock the cassette, just snap apart the lid, turn it round and lock it into the base of the cassette
- They may be opened and closed as many times as required



### Specifications

<b>Colour</b>	Pink
<b>Quantity</b>	2000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)