

Deltalab Green Tissue Embedding Cassettes

Product Code: 554147

Brand: Deltalab

Deltalab's Green Tissue Embedding Cassettes are made of acetal polymer and have integral lids.


Each cassette has a small tab on the lid for easy opening and closing.

Features

- Air vents allowing more efficient filling with paraffin
- Safety closure
- With one-piece integral lid
- To lock the cassette, just snap apart the lid, turn it round and lock it into the base of the cassette
- They may be opened and closed as many times as required

Specifications

Colour	Green
Quantity	2000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com