

## Deltalab White Tissue Embedding Cassettes (Slimsette)

**Product Code:** 555141

**Brand:** Deltalab

Made of polyacetal, Deltalab's White Tissue Embedding Cassettes (Slimsette) provide quick drainage thanks to their 114 vents measuring 5 x 1mm. Their hinged lids can be opened and closed as often as necessary; they may be separated as required. They are supplied assembled, saving time and making storage easier.

These cassettes can be easily handled with one hand and closed safely. Their air vents allow more efficient filling with paraffin. The writing area is slanted at a 45° angle and suitable with most labelling instruments.

### Specifications

<b>Dimensions (mm)</b>	41 x 28.5 x 6
<b>Colour</b>	White
<b>Quantity</b>	1500



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)