


Deltalab Aluminium Screw Caps for Screw Threaded Glass Tubes

Product Code: 617000

Brand: Deltalab

Deltalab's Aluminium Screw Caps for Screw Threaded Glass Tubes are of pharmaceutical quality and fully autoclavable.

These caps have an internal liner made of caoutchouc that are suitable for flat bottom soda glass tubes.



Specifications

Suitable Tube	617145
Type of Cap	Aluminium
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com