

## Deltalab 24ml Flat Bottom Soda Glass Screw Threaded Tubes

**Product Code:** 617145

**Brand:** Deltalab

The high quality of Deltalab's 24ml Flat Bottom Soda Glass Screw Threaded Tubes are reflected in the uniformity of their wall thickness and of their diameter and height dimensions.

They are supplied in small quantities per case for a more convenient use in laboratories.

They also come without caps.



### Specifications

<b>Total Capacity (ml)</b>	24
<b>Screw</b>	15 - 415
<b>Internal Diameter (mm), tube</b>	14.20
<b>External Diameter (mm), tube</b>	17.75
<b>Height (mm)</b>	144.5
<b>Thickness (mm)</b>	1.05
<b>Quantity</b>	252

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)