

## Deltalab 150mm Glass Pasteur Pipettes, Open Mouth

**Product Code:** 701

**Brand:** Deltalab

Made of soda glass, Deltalab's 150mm Glass Pasteur Pipettes, Open Mouth enable users to easily reach into deep vials and can be used for studies in various scientific fields (e.g. bacteriology, immunology).



### Specifications

<b>Dimensions</b>	<b>Mouth Diameter (mm)</b>	6.95 ± 0.15
	<b>Tip Diameter (mm)</b>	1.2 ± 0.15
	<b>Thickness (mm)</b>	0.53 ± 0.03
<b>Total Length (mm)</b>		150
<b>Stem Length (mm)</b>		59
<b>Quantity</b>		1000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)