

Deltalab 230mm Glass Pasteur Pipettes, Open Mouth

Product Code: 702

Brand: Deltalab

Made of soda glass, Deltalab's 230mm Glass Pasteur Pipettes, Open Mouth enable users to easily reach into deep vials and can be used for studies in various scientific fields (e.g. bacteriology, immunology).



Specifications

Dimensions	Mouth Diameter (mm)	6.95 ± 0.15
	Tip Diameter (mm)	1.2 ± 0.15
	Thickness (mm)	0.53 ± 0.03
Total Length (mm)		230
Stem Length (mm)		130
Quantity		1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com