

## Ceti Cold Light Source Illuminator with Metal Housing (150W)

**Product Code:** 7063.1502

**Brand:** Ceti

Ceti's Cold Light Source Illuminator with Metal Housing (150W) features a double swan neck arm.



### Specifications

<b>Main Power</b>	220V AC 50Hz (Euro connector)
<b>Light Intensity</b>	Regulation with electric control 0 - 100%
<b>Colour Temperature</b>	3200°K
<b>Heat Filter Bulb</b>	21V 150W Halogen reflector bulb (500 burning hours)
<b>Cooling System</b>	Axial ventilator
<b>Packing Dimensions (W x D x H)</b>	38 x 26 x 21cm

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)