

Ceti Fiber Optic Illuminator (150W)

Product Code: 7063.1555

Brand: Ceti

The Ceti Fiber Optic Illuminator (150W) is used to provide acute, cool illumination for various optical applications.

With its double swan neck arm, this light source is designed to border with fibre optic light guides to deliver a continuous light supply that can be manoeuvred or conducted in many different ways. The CETI Fiber Optic Illuminator (150W) provides easy-to-use, active illumination making it very beneficial for several applications.



Specifications

Main Power	220V AC 50Hz
Light Intensity	Regulation with electronic control 0 - 100%
Colour Temperature	3200°K
Bulb	150W Halogen reflector bulb (estimated 500 burning hours)
Cooling System	Axial ventilator
Packing Dimensions (W x D x H)	38 x 26 x 21cm
Gross Weight (kg)	3

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com