


Deltalab 230mm Glass Pasteur Pipettes, Closed Tip

Product Code: 712

Brand: Deltalab

Made of natural soda glass, Deltalab's 230mm Glass Pasteur Pipettes, Closed Tip have cotton in the mouth. The closed part can be opened easily by pushing it slightly.



Specifications

Dimensions	Mouth (mm)	7
	Tip (mm)	1.4 to 1.6
	Thickness (mm)	0.50
Total Length (mm)		230
Stem Length (mm)		130
Quantity		1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com