

Deltalab 5ml Round Bottom Borosilicate Glass Screw Threaded Tubes

Product Code: 713100

Brand: Deltalab

The high quality of Deltalab's 5ml Round Bottom Borosilicate Glass Screw Threaded Tubes are reflected in the uniformity of their wall thickness and of their diameter and height dimensions.

They are supplied in small quantities per case for a more convenient use in laboratories.

They also come without caps.



Specifications

Total Capacity (ml)	5
Screw	13 - 415
Internal Diameter (mm), tube	10.65
External Diameter (mm), tube	12.75
Height (mm)	100
Thickness (mm)	1.00
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com