

Deltalab 18ml Round Bottom Borosilicate Glass Screw Threaded Tubes

Product Code: 716150

Brand: Deltalab

The high quality of Deltalab's 18ml Round Bottom Borosilicate Glass Screw Threaded Tubes are reflected in the uniformity of their wall thickness and of their diameter and height dimensions.

They are supplied in small quantities per case for a more convenient use in laboratories.

They also come without caps.

Specifications

Total Capacity (ml)	18
Screw	15 - 415
Internal Diameter (mm), tube	14
External Diameter (mm), tube	16
Height (mm)	150
Thickness (mm)	1.05
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com