

## Deltalab 70mm Microhaematocrit Capillary Tubes (With Heparin)

**Product Code:** 7401

**Brand:** Deltalab

Deltalab's 70mm Microhaematocrit Capillary Tubes (With Heparin) are soda neutral glass tubes with colour-coded printing for easier identification with sodium heparin (red).

They are supplied inside a glass tube with a plastic cap, specifying code, lot, and expiry date.

### Specifications

<b>Description</b>	70mm with heparin
<b>External Diameter</b>	1.5 - 1.6
<b>Volume (µl)</b>	70
<b>Quantity</b>	100

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)