


Deltalab 70mm Microhaematocrit Capillary Tubes (With Heparin)

Product Code: 7401

Brand: Deltalab

Deltalab's 70mm Microhaematocrit Capillary Tubes (With Heparin) are soda neutral glass tubes with colour-coded printing for easier identification with sodium heparin (red).

They are supplied inside a glass tube with a plastic cap, specifying code, lot, and expiry date.



Specifications

Description	70mm with heparin
External Diameter	1.5 - 1.6
Volume (µl)	70
Quantity	100

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com