


Deltalab 70mm Microhaematocrit Capillary Tubes (Without Heparin)

Product Code: 7411

Brand: Deltalab

Deltalab's 70mm Microhaematocrit Capillary Tubes (Without Heparin) are soda neutral glass tubes with colour-coded printing for easier identification without heparin (blue).

They are supplied inside a glass tube with a plastic cap, specifying code, lot, and expiry date.



Specifications

| | |
|--------------------------|----------------------|
| Description | 70mm without heparin |
| External Diameter | 1.5 - 1.6 |
| Volume (µl) | 70 |
| Quantity | 100 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com