

Ceti DIGIT-012 Clinical Hand-Held Analogue Refractometer

Product Code: 8100.0120M

Brand: Ceti

The Ceti DIGIT-012 Clinical Hand-Held Analogue Refractometer is used for measuring proteins in serum and urine specific gravity.

Its main glass prism is fixed with a chemical resistant epoxy.

It also features a focusing eyepiece with a soft rubber eyecup.



Specifications

Scale	0 - 12g/dl (Protein) 1.3330 to 1.3600 (Refractive Index) 1.000 to 1.050 (Specific Gravity)
Precision	0.2g/dl 0.0005 nD 0.002sg
Prism Cover	Plexi
Dimensions (mm)	55 x 80 x 225
Weight (g)	160
Calibration	Zero calibration with distilled water at 20°C

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com