

## Ceti DIGIT-015 Sugar/Oil Hand-Held Analogue Refractometer

**Product Code:** 8100.0150M

**Brand:** Ceti

The Ceti DIGIT-015 Sugar/Oil Hand-Held Analogue Refractometer is a pocket refractometer used for measuring low sugar concentrations (except in fruit juices) and low refractive values (except of cutting oils). It is suitable for general applications in the industrial sector and offers precision and versatility.

This laboratory refractometer features a main glass prism fixed with a chemical resistant epoxy and a focusing eyepiece.



### Specifications

<b>Scale</b>	0 - 15% (Brix)
	1.3330 to 1.3560 (Refractive Index)
<b>Precision</b>	0,2% (Brix)
	0.0005 nD
<b>Prism Cover</b>	Plexi
<b>Dimensions (mm)</b>	55 x 80 x 205
<b>Weight (g)</b>	160
<b>Calibration</b>	Zero calibration with distilled water at 20°C

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)