

## Ceti DIGIT-080 Sugar Hand-Held Analogue Refractometer

**Product Code:** 8100.0800M

**Brand:** Ceti

The Ceti DIGIT-080 Sugar Hand-Held Analogue Refractometer is a lightweight sugar pocket refractometer used for measuring high sugar concentrations (except in juices, honey, jam, syrup, canned food).

This laboratory refractometer features a main and secondary glass prism fixed with a chemical resistant epoxy and focusing eyepiece with soft rubber eyecup.



### Specifications

<b>Scale</b>	One scale
	0 - 80% (Brix)
<b>Precision</b>	0,5% (Brix)
<b>Prism Cover</b>	Metal
<b>Dimensions (mm)</b>	55 x 80 x 225
<b>Weight (g)</b>	180
<b>Calibration</b>	Zero calibration with distilled water at 20°C

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)