

## Ceti DIGIT-090 Sugar Hand-Held Analogue Refractometer

**Product Code:** 8100.0900M

**Brand:** Ceti

The Ceti DIGIT-090 Sugar Hand-Held Analogue Refractometer is used for measuring high sugar concentrations.

The main glass prism of this laboratory refractometer is fixed with a chemically resistant epoxy. It also has an adjustment ring for colour dispersion.



### Specifications

<b>Scale (triple)</b>	0 - 42% (Brix)
	42 - 71% (Brix)
	71 - 90% (Brix)
<b>Precision</b>	0,2% (Brix)
<b>Prism Cover</b>	Plexi
<b>Dimensions (mm)</b>	55 x 125 x 235
<b>Weight (g)</b>	900
<b>Calibration</b>	Zero calibration with distilled water at 20°C

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)