

Ceti DIGIT-2862 Sugar Hand-Held Analogue Refractometer

Product Code: 8100.2862

Brand: Ceti

The Ceti DIGIT-2862 Sugar Hand-Held Analogue Refractometer is a lightweight sugar pocket refractometer used for measuring mid-range concentrations (except in canned food or sugar solutions).

This laboratory refractometer features a main glass prism with a chemical resistant epoxy and a focusing eyepiece with a soft rubber eyecup.



Specifications

Scale	28 - 62% (Brix)
Precision	0,2% (Brix)
Prism Cover	Plexi
Dimensions (mm)	55 x 80 x 205
Weight (g)	180
Calibration	Saturated salt solution at 29.6%

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com