

## Ceti DIGIT-30-I40 Oechsle Hand-Held Analogue Refractometer

**Product Code:** 8100.3014

**Brand:** Ceti

The Ceti DIGIT-30-I40 Oechsle Hand-Held Analogue Refractometer is used for measuring sugar or alcohol volumes in non-fermented fruit or during wine production.

This refractometer has a main glass prism fixed with a chemically resistant epoxy and a focusing eyepiece with a soft rubber eyecup.



### Specifications

<b>Scale</b>	0 - 32% (Brix)
	30 - 140% (Oechsle)
<b>Precision</b>	0.20%
	0.10%
<b>Prism Cover</b>	Plexi
<b>Dimensions (mm)</b>	55 x 80 x 205
<b>Weight (g)</b>	125
<b>Calibration</b>	Zero calibration with distilled water at 20°C

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)