

Ceti DIGIT-5890 Honey Hand-Held Analogue Refractometer

Product Code: 8100.5890M

Brand: Ceti

The Ceti DIGIT-5890 Honey Hand-Held Analogue Refractometer is used for measuring high concentrations (ex. in honey, jam, marmalade).

This refractometer has a main glass prism fixed with a chemically resistant epoxy and a focusing eyepiece with a hard plastic eyecup.



Specifications

Scale	58 - 90% (Brix)
	12 - 27% (Water)
	38 - 43% (Baumé)
Precision	0,5% (Brix)
	1% (Water)
	0,5% (Baumé)
Prism Cover	Plexi
Dimensions (mm)	55 x 80 x 205
Weight (g)	250
Calibration	Calibration fluid

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com