

Ceti DIGIT-010 ATC Sugar Hand-Held Analogue Refractometer

Product Code: 8145.0001

Brand: Ceti

The Ceti DIGIT-010 Sugar Hand-Held Analogue Refractometer is a sugar pocket refractometer, which measures low concentrations of sugar or cutting oil(s) with automatic temperature compensation.



Specifications

Scale	0 - 10% (Brix)
Precision	0,1% (Brix)
Prism Cover	Plexi
Dimensions (mm)	55 x 80 x 225
Weight (g)	140
Calibration	Zero calibration with distilled water between 10 and 30°C

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com