

## Ceti DIGIT-100 ATC Salinity Hand-Held Analogue Refractometer

**Product Code:** 8145.0100M

**Brand:** Ceti

The Ceti DIGIT-100 ATC Salinity Hand-Held Analogue Refractometer is a light and simple salinity pocket refractometer used for measuring low salinity concentrations of natural or artificial salt in water with automatic temperature compensation.

This laboratory refractometer features a main glass prism fixed with a chemical resistant proxy and a focusing eyepiece.



### Specifications

<b>Scale</b>	0 - 100% (Salinity)
	1.000 to 1.070sg (Gravity)
<b>Precision</b>	1% (Salinity)
	0.001sg
<b>Prism Cover</b>	Plexi
<b>Dimensions (mm)</b>	55 x 80 x 225
<b>Weight (g)</b>	250
<b>Calibration</b>	Zero calibration with distilled water between 10 and 30°C

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)