

Ceti DIGIT-025 ATC Alcohol Hand-Held Analogue Refractometer

Product Code: 8145.0250M

Brand: Ceti

The Ceti DIGIT-025 ATC Alcohol Hand-Held Analogue Refractometer is used for measuring alcohol volume in non-fermented fruits and during wine production or vintages.

This refractometer has a main glass prism fixed with a chemically resistant epoxy and a focusing eyepiece with a soft rubber eyecup.



Specifications

Scale	0 - 25% (Alcohol)
Precision	0.2%
Prism Cover	Plexi
Dimensions (mm)	55 x 70 x 200
Weight (g)	250
Calibration	Zero calibration with distilled water at 20°C

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com