

Ceti DIGIT-032 ATC Sugar Hand-Held Analogue Refractometer

Product Code: 8145.0320M

Brand: Ceti

The Ceti DIGIT-032 ATC Sugar Hand-Held Analogue Refractometer is a portable sugar pocket refractometer used for measuring low sugar concentrations (except in fruit and soft drinks) with automatic temperature compensation.

This laboratory refractometer features a main glass prism fixed with a chemical resistant proxy and a focusing eyepiece.



Specifications

Scale	0 - 25% (Alcohol/Volume)
Precision	0,2% (Brix)
Prism Cover	Plexi
Dimensions (mm)	55 x 80 x 205
Weight (g)	210
Calibration	Zero calibration with distilled water between 10 and 30°C

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com