

## Ceti DIGIT-5890 ATC Honey Hand-Held Analogue Refractometer

**Product Code:** 8145.5890M

**Brand:** Ceti

The Ceti DIGIT-5890 ATC Honey Hand-Held Analogue Refractometer is a compact and durable honey pocket refractometer that provides simple and accurate Brix readings in honey.

This lightweight laboratory refractometer features a main glass prism fixed with a chemical resistant proxy and a focusing eyepiece with soft rubber eyecup.



### Specifications

<b>Scale</b>	58 - 90% (Brix)
	12 - 27% (Water)
	38 - 43% (Baumé)
<b>Precision</b>	0,5% (Brix)
	1% (Water)
	0,5% (Baumé)
<b>Prism Cover</b>	Plexi
<b>Dimensions (mm)</b>	55 x 80 x 205
<b>Weight (g)</b>	250
<b>Calibration</b>	Zero calibration with calibrated fluid (50% Brix between 10 and 30°C)

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)