


Deltalab 27ml Round Bottom Soda Glass Tubes

Product Code: 816160

Brand: Deltalab

The high quality of Deltalab's 27ml Round Bottom Soda Glass Tubes are reflected in the uniformity of their wall thickness and of their diameter and height dimensions.

They are supplied in small quantities per case for a more convenient use in laboratories.



Specifications

Total Capacity (ml)	27
Internal Diameter, tube (mm)	14.40
External Diameter, tube (mm)	16
Height (mm)	160
Thickness (mm)	0.55
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com