


## Deltalab 28ml Round Bottom Soda Glass Tubes

**Product Code:** 818150

**Brand:** Deltalab

The high quality of Deltalab's 28ml Round Bottom Soda Glass Tubes are reflected in the uniformity of their wall thickness and of their diameter and height dimensions.

They are supplied in small quantities per case for a more convenient use in laboratories.



### Specifications

|                                     |      |
|-------------------------------------|------|
| <b>Total Capacity (ml)</b>          | 28   |
| <b>Internal Diameter, tube (mm)</b> | 15   |
| <b>External Diameter, tube (mm)</b> | 18   |
| <b>Height (mm)</b>                  | 150  |
| <b>Thickness (mm)</b>               | 0.85 |
| <b>Quantity</b>                     | 500  |

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)